Datasheet

 Mouse mAb to
 Keratin 4,5,6,8,10,13,18

 Clone
 C-11

 Isotype
 IgG1-к

Source

A BALB/c mouse was immunized with a cytoskeleton preparation of human A431 carcinoma cell line. Fusion partner: P3-X63-Ag8.653.

Specifications

C11 reacts with keratins : 4,5,6,8,10,13 and 18. This is a broad-spectrum antibody, which has been reported to differentiate epithelial tumors from non-epithelial tumors. Many studies have shown the usefulness of keratins as markers in cancer research and tumor diagnosis.



Positive: cat, chicken, cow, dog, ferret, goat, guinea pig, human, kangaroo, marmoset, monkey, mouse, rabbit, rat, *Xenopus laevis*.



Figure 1: Colon carcinoma stained with C11 (paraffin)

Applications

Immunocytochemistry on frozen and paraffin sections. Immunofluorescence tests. Strong reaction in immunoblotting.

Flow cytometry	Frozen sections	Immunofluorescence	Paraffin sections	Western blot
+	+	+	Citrate/protease	+

Format

Produced in tissue culture, contains no host Ig. Antibodies are affinity purified and presented in PBS with 0,02 % sodium azide.

Stored at 4°C- 8°C, shelf life is at least 24 months after purchase.

Dilution advice

- Flow Cytometry (0.5-1.0 μg/million cells in 0.1 ml, fix cells in 4% PFA for 10 min, at 4°C, permeabilize with 0,2% saponin or digitonin for 15 min, at 4°C).
- > Immunoblotting (1 μ g/ml for 2h at RT).
- Immunofluorescence (0.5-1 μ g/ml).
- Immunohistology (formalin-fixed: 1-2 μg/ml for 30 min at RT; requires boiling tissue sections in 10mM citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes or digestion of sections with trypsin at 1mg/ml PBS, 15 min at RT).

Positive control

A431 cells, skin, lung.

Datasheet



References

- Vojtěšek, B. et al. *Folia Biol.*, **35(6)**: 373-382, (1989)
 Vojtěšek, B. et al. *Neoplasma* **37(3)**: 333-342 (1990).
- Kovařík J. et al. Int. J. Cancer, Suppl. 3: 50-55, (1988).
- Kovařík J. et al. J. Tumor Marker Oncol. 5: 219 (1990).