Datasheet

Blood Group A
(CD173)
HE-10
IgM-к

Source

A BALB/c mouse was immunized with a mixture of erythrocytes of blood group A1 and a glycoprotein fraction isolated from the saliva of secretors with blood group A of human origin. Fusion partner: P3-X63-Ag8.653.

Specifications

HE-10 preferably reacts with determinants of chain A type 3 and 4 and chain H type 3 and 4, but not with type 1 and 2 chain structures. It is not reactive with immunodominant A trisaccharide. Increased expression of this antigen has been observed on some tumor tissues such as gastric carcinomas, urothelial carcinomas, and colon carcinomas. HE-10 does not react with normal tissue sections of donors with blood group B and 0 but it reacts specifically with malignant tissues in these individuals.



Figure 1: Human rectum stained with HE-10 (paraffin)

Species reactivity

Positive: human.

Applications

HE-10 is applicable for red cell agglutination, tissue staining and immunofluorescence tests.

Agglutination	Flow cytometry	Frozen sections	Immunofluorescence	Paraffin sections
+	+	+	+	Citrate

Format

Produced in tissue culture, contains no host Ig. Antibodies are affinity purified and presented in PBS with 0,02% sodium azide.

Stored at 4°C- 8°C, shelf life is at least 24 months after purchase.

Dilution advice

- > Agglutination
- Flow cytometry (0,5–1,0 μg/million cells in 0,1 ml).
- > Immunofluorescence (0,5–1 μ g/ml).
- Immunohistology (formalin-fixed: 1-2 μg/ml for 30 min. at RT; requires boiling tissue sections in 10mM citrate buffer, pH 6,0, for 10-20 min followed by cooling at RT for 20 min).

Positive control

KG1 cells or human colorectal carcinoma.

Datasheet



References

- \succ
- ≻
- Němec M. et al. *Vox Sang* **52**:125-8 (1987). Vanák J. et al. *Neoplasma* **36**(**4**): 479-487 (1989). Tichý, M. et al. *Neoplasma* **37**(**4**): 451-459 (1990). \triangleright
- > Tichý, M. et al. Acta Univ Palacki Olomuc Fac Med. 126: 57-69 (1990).